OSMOCELL

CE 0373

## ANTICELLULITE

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CE 0373

## ANTICELLULITE

Cellulite, which is more precisely defined as oedematous fibrosclerotic panniculopathy (E.F.P.), results from adermal and hypodermal alteration. In optimal conditions, the adipose cells (adipocytes) that are normally present in subcutaneous tissue act as an energy reserve for the body, which burns fat whenever it needs fuel. In case of stasis of the venous circulation, this “reserve” becomes hard to use and builds up to the point of compressing capillaries which, fragile by nature, start transuding plasma through their walls, which have become porous. Plasma infiltrates between the cells and, in the course of time, causes adipose tissue inflammation with the formation of subcutaneous tissue fibrosis. Capillaries are further compressed, their walls become more porous and permeable, allowing the outflow of water, which invades interstitial space between adipocytes (adipose cells). In the course of time, drainage of excess water becomes increasingly difficult and triggers a “vicious circle” that further aggravates the disease.


Medical Device  
Class IIa

## SOLUTION

The objective was to intervene on the principal characterising factors of cellulite, namely fluid stasis and structural alteration of adipose tissue. Hence, a strongly hypertonic emulsion was developed with an appropriate partition coefficient between the hydrophilic phase and the lipophilic one.

Said coefficient allows sodium chloride, whose concentration is higher (13%), to permeate the horny layer and recall, by osmosis, the water that is trapped in interstitial spaces between adipose cells, thus encouraging its subsequent elimination through the microcirculation and reducing the inflammation of adipose tissue itself. The presence of caffeine in our formula also ensures an effective lipolytic action that can reduce the adipose ridges that are already present.


OSMOCELL

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The Department of Pharmaceutical Sciences of the University of Milan was requested to conduct a study to evaluate the actual capacity of the emulsion to encourage elimination of water through the skin. The purpose of the study was to evaluate “ex-vivo” the quantity of water absorbed by the hypertonic emulsion vs. a placebo emulsion.

### RESULT

The results prove that, unlike the placebo emulsion, the hypertonic emulsion successfully moves the water placed under the membrane outwards towards the skin surface. One hour after application on the membrane in the adopted conditions, the water content of the hypertonic emulsion increases by 5%, while that of the placebo emulsion drops by about 4%.

EMC italia

Clinical study performed on 30 patients with Edematous Fibro Sclerotic Panniculopathy, with measurements after 28 and 56 days of treatment.

MEASUREMENT	T=28 DAYS	T=56 DAYS
thigh circumference	- 1,44 cm	- 1,88 cm
Panel thickness adipose (ultrasound)	-3,60 mm	-4,20 mm
Skin elasticity (Elastometry)	+ 1,85%	+ 2,53%